

LLLT Series

Wearable neck laser therapy device

650nm

Wearable neck laser therapy device is convenient and easy to use. Red light relieves pain with gentle massage through laser diodes. Compared with other methods, it leads energy to the uncomfortable areas directly. Free to change continuous on and pulse mode. Spend 30 min. a time for promoting circulation and relieve discomfort. Use it continuously for daily healthcare.



Feature

Relaxation and Comfort, Warm Massage

Developed to relieve soreness

- Relieve discomfort
- Alleviates fatigue
- Balance Body and Mind

I Specification

Wavelength	650 ±20nm
Product size	L190/W160/H53 ±5mm
Product Weight	348 ±10g
Battery Capacity	Lithium polymer battery 1900mAh
Laser Safety Classification	Class 1
Maximum Optical Power	0.77mW
Pulse Duration	1 second
Pulse Repetition Rate	1Hz
Operating mode	continuous wave (CW)/ pulse wave (PW)
Laser standard	IEC 60825-1:2014
Operating environment	Temperature 5~40 degrees; Humidity 0~80%
Transportation and storage environment	Temperature 5~40 degrees; Humidity 0~80%

【 Cautions 】

- Avoid light exposure to your eyes directly for protect your eyes.
- Transformer is prohibited for use which may cause dangerous event.
- Do not use fast charging for this product, as it may cause danger.
- Please keep clean on the area of laser output, prevent the dirties affect energy transmission when use.
- Please keep it away from the touch by children or pets in order to avoid eating or touch components accidentally.
- Do not use on pregnant women, who with skin disease, malignant tumor, wear pacemaker, initial wound within 72 hours, physical abnormalities and abnormal sensory nerves. Avoid to use it during the skin is sensitive, fragile and in rash or wounds condition. Please stop using it immediately and consult doctor if you feel any discomfort.
- Don' t use it with electronic medical equipment in order to avoid the malfunction or dangerous accidents.
- Please store it in dry, cool place and away from water. Please pay attention on moisture and falling during use.
- Remove cloth between device and skin in order to avoid reducing the light intensity.

I LECC TECHNOLOGY CO., LTD.

2F No.2, Ziqiang 3rd Rd., Zhongli Dist., Taoyuan City, 320023, Taiwan.
 TEL | 03 4511967 VAT Number | 16708380
 FAX | 03 4511923 Email | sales@lecc.com.tw

I SHANGHAI LECC OPTO CO., LTD.

No.1939, Da Ye Rd., Feng-Xian Dist., Shanghai City, 201402, China.
 TEL | 021 57407400 TAX Number | 913100007630197034
 FAX | 021 57407403 Email | sales@lecc.com.tw