

LLLT Series

Nasal laser therapy device

650nm

Nasal laser therapy device promotes nasal cavities.
Improves the energy of the nasal cavities and relieves dark circles.
Improves circulation, say goodbye to dark circles.



Feature

- Olfactory conditioning
- Relaxation of nasal muscles
- Compact and portable
- Convenient and safe to use anytime, anywhere

I Specification

Wavelength	650 ±20nm
Product size	L50/W15 ±3mm(Doesn't include wire length)
Product Weight	20 ±1g
Working Voltage	3-5V, 40mA(MAX.)
Laser Safety Classification	Class 2M
Maximum Optical Power	<1mW
Operating mode	continuous wave (CW)
Laser standard	IEC 60825-1:2014
Operating environment	Temperature 10~40 degrees; Humidity 30~75%
Transportation and storage environment	Temperature -20~40 degrees; Humidity 10~90%

【 Cautions 】

- Avoid light exposure to your eyes directly for protect your eyes.
- Transformer is prohibited for use which may cause dangerous event.
- Do not use fast charging for this product, as it may cause danger.
- Please keep clean on the point of laser output and nasal cavities, prevent the dirties affect energy transmission when use.
- Please keep it away from the touch by children or pets in order to avoid eating or touch components accidently.
- Do not use on pregnant women, who with skin disease, malignant tumor, wear pacemaker, initial wound within 72 hours, physical abnormalities and abnormal sensory nerves. Avoid to use it during the skin is sensitive, fragile and in rash or wounds condition. Please stop using it immediately and consult doctor if you feel any discomfort.
- Don' t use it with electronic medical equipment in order to avoid the malfunction or dangerous accidents.
- Please store it in dry, cool place and away from water. Please pay attention on moisture and falling during use.

I LECC TECHNOLOGY CO., LTD.

2F No.2, Ziqiang 3rd Rd., Zhongli Dist., Taoyuan City, 320023, Taiwan.
 TEL | 03 4511967 VAT Number | 16708380
 FAX | 03 4511923 Email | sales@lecc.com.tw

I SHANGHAI LECC OPTO CO., LTD.

No.1939, Da Ye Rd., Feng-Xian Dist., Shanghai City, 201402, China.
 TEL | 021 57407400 TAX Number | 913100007630197034
 FAX | 021 57407403 Email | sales@lecc.com.tw